



Coffalyser digitalMLPA™ Release Notes



General information

Coffalyser digitalMLPA is Research Use Only (RUO) software designed for the analysis of digitalMLPA data generated using a SALSA® digitalMLPA probemix* as described in the digitalMLPA General Protocol

* To be used in combination with a SALSA® digitalMLPA Reagent Kit and SALSA® Barcode Solutions Plates.

Version numbers

Software version 221020.1234

Reporting

Module versions:

 Configuration
 2022-10-20 12:34

 Data conversion
 2021-10-19 09:58

 Fragment analysis
 2021-10-19 09:58

 Comparative analysis
 2021-10-19 09:58

 Extended analysis
 n/a

www.mrcholland.com page 1 of 6

2021-10-19 09:58



Changes in Coffalyser digitalMLPA compared to previous versions

Table 1 Overview of changes per release

		Summary of changes		
Release	v221020.1234	Product sheets for lots D006-A1-1022 and D007-A1-0922 added to configuration module		
	v220221.1522	Product sheet for lot D001-C1-1221 added to configuration module		
	v211020.1245	Changes in data processing, quality checking and reporting as well as implementation of bug fixes to improve performance. Product sheet for lot D001-B1-1119 updated.		
	v210729.1257	Product sheet for lot D001-B1-0421 and a restricted product sheet added to configuration mo		
	v201211.1315	Restricted product sheet added to configuration module		
	v201007.1443	Restricted product sheet added to configuration module		
	v200609.0921	First external release		

Table 2 Overview of changes per probemix lot per release

		Product sheets						
		D001				D007		
		B1-1119	B1-0421	C1-1221	A1-1022	A1-0922		
Release	v221020.1234	No changes	No changes	No changes	New	New		
	v220221.1522	No changes	No changes	New	-	-		
	v211020.1245	Data processing, quality checking and reporting updated. Updated product sheet.	Data processing, quality checking and reporting updated. Product sheet unchanged.	-	-	-		
	v210729.1257	No changes	New	-	-	-		
	v201211.1315	No changes	-	-	-	-		
	v201007.1443	No changes	-	-	-	-		
	v200609.0921	New	-	-	-	-		

We recommend to always use the latest version of Coffalyser digitalMLPA. For software and product sheet updates in bold, (re)validation is recommended.

www.mrcholland.com page 2 of 6



Software version 221020.1234

Configuration: 2022-10-20 12:34 | Data conversion: 2021-10-19 09:58 | Fragment analysis: 2021-10-19 09:58 | Comparative analysis: 2021-10-19 09:58 | Extended analysis: n/a | Reporting: 2021-10-19 09:58.

General

- Configuration module version modified to accommodate new commercially available digitalMLPA probemixes.
- Product sheets for previous limited-release probemix versions removed from the configuration module.

Probemix-specific

D006-A1:

- Product sheet released for lot D006-A1-1022 for research use only (RUO).

D007-A1:

Product sheets released for lot D007-A1-0922 for research use only (RUO): D007-A1 default product sheet and D007-A1 near haploid_low hypodiploid product sheet.

D001-C1:

Product sheet for lot D001-C1-1221 remains unchanged.

D001-B1:

- Product sheets for lot D001-B1-1119 and D001-B1-0421 remain unchanged.

Software version 220221.1522

Configuration: 2022-02-21 15:22 | Data conversion: 2021-10-19 09:58 | Fragment analysis: 2021-10-19 09:58 | Comparative analysis: 2021-10-19 09:58 | Extended analysis: n/a | Reporting: 2021-10-19 09:58.

General

Configuration module version modified to accommodate commercially available digitalMLPA probemixes.

Probemix-specific

D001-C1:

- Product sheet released for lot D001-C1-1221 for research use only (RUO).

D001-B1:

- Product sheets for lot D001-B1-1119 and D001-B1-0421 remain unchanged.

Software version 211020.1245

Configuration: 2021-10-20 12:45 | Data conversion: 2021-10-19 09:58 | Fragment analysis: 2021-10-19 09:58 | Comparative analysis: 2021-10-19 09:58 | Extended analysis: n/a | Reporting: 2021-10-19 09:58.

General

Data processing:

- An option is included to ignore undefined samples (barcodes) present in the FASTQ file(s) and exclude these from analysis.
- Assignment of reads to probes is updated to avoid the possibility of reads being assigned to more than one probe, which would prevent the analysis from continuing.

www.mrcholland.com page 3 of 6



- Coffa files generated by previous releases are not compatible with this release. Any re-analysis should start from the FASTQ file(s).
- Inter-normalisation of reference samples is updated.
- The number of reads required per sample when the sample was not included in the experiment definition has been increased to avoid the creation of unnecessary undefined experiment(s).
- Barcode plates 3 and 4 are added to allow analysis of experimental data in which SALSA® Barcode Solutions Plates 3 and 4 are used.

Quality checking:

- A warning is added to the Quality table in the report in case less than three reference samples are included in an experiment with defined reference samples.
- A warning is included to the Quality table in the report in case there are probes present in the sample that cannot be classified due to inadequate normalisation.
- Reference Sample Quality problems will be penalised more severely, and issues will therefore more easily result in sample quality failure.
- DNA input estimation calculation updated.
- Processing of depurination sensitive control probes is updated.

Reporting:

- The Overview.pdf report has been updated:
 - Various small textual and lay-out changes.
- The General Details.pdf report has been updated:
 - Product selection field updated and renamed to Product detection.
 - Quality Table updated:
 - New checks on product detection, probes not being classified and number of reference samples added.
 - Read depth (qualified) and disqualified lot identification probes removed.
 - Order of checks updated.
 - Various small textual and lay-out changes.
- The Sample Results.pdf report has been updated:
 - o Product selection field updated and renamed to Product detection.
 - O Quality Table updated:
 - New checks on product detection, probes not being classified and number of reference samples added.
 - Read depth (qualified) and disqualified lot identification probes removed.
 - Order of checks updated.
 - o Colouring of *Inter ratio*, *Within normal range* and *Different from reference population* columns changed when data is incomplete.
 - o Various small textual and lay-out changes.
- The No Dna.pdf report has been updated:
 - Various small textual and lay-out changes.
- The Ratios.xlsx report has been updated:
 - o File format changed from .xls to .xlsx.
 - o Additional tab containing Reference Sample Quality data added.
 - Ratios in the Inter ratios sheet with missing classification are highlighted in grey.
 - o Various small textual and lay-out changes.

Bug fixes:

- Version number of Configuration as displayed on reports is made independent of local time settings of computer.
- Memory usage is reduced during report generation to prevent out-of-memory error.
- Regional reporting issues are fixed.
- In Command Line Interface additional checks are included to ensure the same Product Sheet is used during the complete analysis.

Other:

- Required .NET Framework version is changed from 4.7 to 4.8.
- Link to Coffalyser digitalMLPA User Manual on website is added to the start screen.

www.mrcholland.com page 4 of 6



Various small textual and lay-out changes in the Graphical User Interface.

Related to the Coffalyser digitalMLPA Release Notes:

- Intended purpose on first page of these Release Notes renamed to General information as the Coffalyser digitalMLPA software is intended for Research Use Only.
- SALSA® Barcode Solutions updated to SALSA® Barcode Solutions Plates.

Probemix-specific

D001-B1:

- Product sheet for lot D001-B1-0421 remains unchanged.
- Product sheet for lot D001-B1-1119 contains updated information.

Software version 210729.1257

Configuration: 2021-07-29 12:57 UTC | Data conversion: 2020-05-12 20:10 | Fragment analysis: 2020-05-12 20:10 | Comparative analysis: 2020-05-12 20:10 | Extended analysis: n/a | Reporting: 2020-05-12 20:10.

General

 Configuration module version modified to accommodate commercially available and limitedrelease[#] digitalMLPA probemixes.

Probemix-specific

D001-B1:

- Product sheet released for lot D001-B1-0421 for research use only (RUO).
- Product sheet for lot D001-B1-1119 remains unchanged.

Software version 201211.1315

Configuration: 2020-12-11 13:15 UTC | Data conversion: 2020-05-12 20:10 | Fragment analysis: 2020-05-12 20:10 | Comparative analysis: 2020-05-12 20:10 | Extended analysis: n/a | Reporting: 2020-05-12 20:10.

General

- Configuration module version modified to accommodate a limited-released# digitalMLPA probemix. This change does not affect analysis of the D001-B1 probemix.

Software version 201007.1443

Configuration: 2020-10-07 14:43 UTC | Data conversion: 2020-05-12 20:10 | Fragment analysis: 2020-05-12 20:10 | Comparative analysis: 2020-05-12 20:10 | Extended analysis: n/a | Reporting: 2020-05-12 20:10.

General

Configuration module version modified to accommodate a limited-released[#] digitalMLPA probemix. This change does not affect analysis of the D001-B1 probemix.

Software version 200609.0921

Configuration: 2020-06-09 07:21 UTC | Data conversion: 2020-05-12 20:10 | Fragment analysis: 2020-05-12 20:10 | Comparative analysis: 2020-05-12 20:10 | Extended analysis: n/a | Reporting: 2020-05-12 20:10.

General

External release.

www.mrcholland.com page 5 of 6



Probemix-specific #

D001-B1:

- Product sheet released for research use only (RUO).

www.mrcholland.com page 6 of 6

[#] For details about changes made to analysis sheets that are available for limited-release digitalMLPA probemixes, please contact <u>info@mrcholland.com</u>.